



# Arisure<sup>®</sup> Closed System Transfer Device

FDA Cleared Under ONB Product Code



Keeps you safe so you can focus on what matters – the **PATIENT**.



## Arisure<sup>®</sup> CSTD is Dry by Design\*

A needle-free, dual membrane system to enable airtight and leakproof transfer of hazardous medications.

**SAFE | SIMPLE | SECURE**

# ARISURE® CLOSED VIAL ADAPTER (CVA)

The Arisure Closed Vial Adapter (CVA) is a safe and easy-to-use vial access device for use with hazardous medications. The check valve design allows for passive pressure equalization eliminating the need to prime the vial.

## SAFE

Needle-free, dual membrane, airtight and leakproof technology. Mechanically closed to prevent the escape of harmful vapor and aerosols into the environment, and microbiologically closed maintaining drug vial sterility for up to 7 days.

## SIMPLE

Low attachment force and non-coring spike help eliminate problems during vial access.

## SECURE

Locking technology to ensure secure vial to device connection.

## PRODUCT DETAILS

- 13mm, 20mm, and 28mm variations available to fit all vial sizes
- Needle-free technology
- Vapor capture expansion chamber



**SAFE: MADE TO PROTECT YOU & YOUR TEAM**

# ARISURE® CLOSED MALE LUER (CML)

The Arisure CML is easy-to-use and completely needle-free.

## SAFE

Airtight and leakproof technology to protect healthcare workers and patients from hazardous drug administration.

## SIMPLE

Can easily integrate into current practice, and help to ensure compliance. Compact design for ease-of-use, with fewer parts, pieces, and steps.

## SECURE

Ratcheting technology, CML locks into place creating a passive system to enhance safety for nursing ensuring they are protected.

## PRODUCT DETAILS

- Forms an airtight and leakproof connection with Arisure Neutral Valve (ANV)
- Remains closed when not accessed



## ARISURE® SYSTEM ADVANTAGE - DRY BY DESIGN\*

The Arisure CML is compatible with \*\*BD™ SmartSite™ and \*\*\*Baxter ONE-LINK & CLEARLINK connectors to achieve a closed system with needle-free, dual membrane, airtight and leakproof technology. Integrate Arisure as a passive system into your current practice without the need for additional SKUs and special connectors with fewer steps to streamline to USP 800 compliance.

# ARISURE® NEUTRAL VALVE (ANV)

The Arisure Neutral Valve (ANV) is a safe and easy-to-use needle-free access valve used for the transfer of potentially hazardous fluids and medications using an intuitive luer lock.

*\*Data on file with manufacturer.*

*\*\*SmartSite is a trademark of Becton, Dickinson and Company, its subsidiaries or affiliates. Yukon Medical Inc. is not affiliated with Becton, Dickinson and Company.*

*\*\*\*ONE-LINK AND CLEARLINK are trademarks of Baxter International Inc, its subsidiaries or affiliates. Yukon Medical Inc. is not affiliated with Baxter International Inc.*



## SIMPLE: MADE FOR MINIMAL WORKFLOW DISRUPTION

## ARISURE® DRY SPIKE (DS)

The Arisure Dry Spike (DS) Non-Vented provides easy, closed method transfer of hazardous medications into a standard intravenous (IV) fluid bag or container.

### SAFE

Needle-free, dual membrane, airtight and leakproof technology that reduces exposure to harmful drug residue when used with the Arisure Closed Male Luer (CML).

### SIMPLE

Adapts into current processes to optimize workflow.

### SECURE

The double seal design allows for secure and safe connections with a variety of IV set spikes.

### PRODUCT DETAILS








- Closed system prevents accidental spills
- Needle-free access port
- Pierceable septum
- Two dedicated lumens within the spike for high flow rate



**SECURE: MADE TO MINIMIZE EXPOSURE TO VAPORS**

The Arisure® Closed System Transfer Device was made for you. It is designed to protect you and your team from exposure to harmful vapors and residue when compounding, transporting, and delivering hazardous drugs.

**The Arisure CSTD is completely needle-free and is appropriate for use in both pharmacy and nursing.**

Image	Item #	Product Name	Priming Volume	Box Qty
	YM101	13mm Arisure CVA with Lubricated Spike	0.11 mL	50
	YM102	20mm Arisure CVA - High Capacity (up to 100 mL)	0.11 mL	50
	YM054	28mm Arisure CVA - High Capacity (up to 100 mL)	0.12 mL	50
	YM104	Arisure Closed Male Luer with Cap	~0.05 mL	100
	YM004	Arisure Neutral Valve (Female Luer)	0.13 mL	100
	YM050	Arisure Vented Dry Spike	<1.5 mL	50
	YM064	Arisure Non-Vented Dry Spike	<1.5 mL	50

## DISTRIBUTION AND ORDERING



YMIV0001revAUG23

**Exclusive U.S. Distributor:**

Progressive Medical, Inc. (PMI)  
997 Horan Drive • Fenton, MO 63026

**Ordering Information:**

Tel: 314.961.5786 • Fax: 314.961.4535  
cs@progressivemedinc.com • [progressivemedinc.com](http://progressivemedinc.com)



Yukon Medical is a leading developer of innovative pharmaceutical preparation and drug delivery devices.

Progressive Medical, Inc. (PMI) is a market leading National Specialty Medical Device Distributor.



Visit us to Learn More about the  
Arisure<sup>®</sup> Closed System Transfer Device

[yukonmedical.com](http://yukonmedical.com)