



# VERSATILE NEEDLE-FREE MEDICATION TRANSFER DEVICE



West Vial Adapter™ transfer devices can be an efficient solution for transfer of medications between vials and syringes.

This versatile device is an efficient needle-free alternative for delivering medications to the sterile field as recommended by AORN<sup>1,2</sup>.

1 AORN: Association of periOperative Registered Nurses; Guidelines for Perioperative Practice: Medication Safety (2019)  
2 Applicable to the perioperative setting only

## Vial Adapter™ transfer devices

- Provides aseptic, needle-free drug preparation
- Available in 13mm and 20mm sizes
- Fewer steps, parts and sharps
- Easy to use with siliconized spike to reduce insertion force
- Reduces end user variability associated with traditional needle aspirations

## Distribution and Ordering

**DISTRIBUTED BY: Progressive Medical, Inc.**

997 Horan Dr. • Fenton, MO 63026

### ORDERING INFORMATION:

Phone: 314-961-5786 • Fax: 314-961-4535

cs@progressivemedinc.com



## Ordering Information

Item#	Description	Quantity
36098131	13mm - Sterile Drug Transfer Device	150/box
36098090	20mm - Sterile Drug Transfer Device	150/box

### Rx Only

The Vial Adapter™ transfer device is a single-use, sterile, medical device intended for the transfer of drugs contained in vials.

Vial Adapter™, Vented Vial Adapter™, Swabable Vial Adapter™, and Swabable Vented Vial Adapter™ transfer devices carry the CE mark (0344). Vial Adapter™, Vented Vial Adapter™, and Swabable Vial Adapter™ transfer devices are 510(k) cleared by the United States Food and Drug Administration. Products are shown for INFORMATION purposes only and may not be approved for marketing in specific regions. Distribution and use are subject to applicable regulatory approvals and requirements for medical devices. The Vial Adapter™ and Vented Vial Adapter™ transfer devices are configurable and may not be suitable for use with all drugs. The Swabable Vial Adapter™ transfer device and the Swabable Vented Vial Adapter™ transfer device inhibit microbial ingress for up to 24 hours and 72 hours, respectively, when used in accordance with the product instructions for use and proper aseptic technique. The use of the Swabable Vial Adapter™ and the Swabable Vented Vial Adapter™ transfer devices should not be interpreted as modifying, extending, or superseding a drug manufacturer's labeling recommendations for storage and expiration dating. The Swabable Vial Adapter™ and the Swabable Vented Vial Adapter™ transfer devices are not compatible for use with all drug products. Refer to drug manufacturer's labeling and use instructions for device configuration compatibility and for shelf life and sterility information of opened or reconstituted drug product. Failure to follow product instructions for use may result in compromised sterility; contamination; leakage (including possible exposure to medication); and/or inadequate medication reconstitution, transfer, and/or dosing. Product misuse could potentially lead to user and/or patient infection, suboptimal therapy, and/or delay in therapy. Please contact your West Pharmaceutical Services, Inc. (West) representative for product availability. Important product and safety information and warnings at: <https://www.westpharma.com/products/vial-adapter-systems/vial-adapters>.

## FOR MORE INFORMATION VISIT WESTPHARMA.COM

West and the diamond logo are registered trademarks of West Pharmaceutical Services, Inc. in the United States and other jurisdictions. Vial Adapter is a trademark of West Pharma. Services IL, Ltd., a subsidiary of West Pharmaceutical Services, Inc.